

**United
Scientific®**

Certificate of Quality

United Scientific certifies that this product has been manufactured in compliance with all applicable procedures and specifications.

Item #: UNST102

Product Description: SteriTrue™ Sterile Sampling Scoops, PS, 250 mL (8 oz) 100 Pk/100

Component Material: Polystyrene (PS)

Expiration Date: Available upon request

Lot Number: Available upon request

Raw Material:

- Heavy Metals Free
- BPA and Phthalate Free
- BSE/TSE compliant

Sterilization:

Products have been sterilized by Gamma radiation.

FDA 21CFR Food Safety:

This product meets the FDA criteria in 21 CFR 177.1520(2017) for contact applications.

RoHS:

Restriction of Hazardous Substance is Electric and Electronic Equipment (RoHS) : Based on the available documentation from raw materials suppliers, this product complies with the requirements of the directives 2011/65/EU, as amended, and 2015/863/EU concerning the limits of Lead (Pb), Mercury (Hg), Hexavalent Chromium (Cr6+), Polybrominated biphenyl(s) (PBB), polybrominated diphenyl ether(s) (PBDE), Bis(2-ethylhexyl) phthalate (DEHP) , Butyl benzyl phthalate (BBP) , Dibutyl phthalate (DBP), Isobutyl phthalate (DIBP) , cadmium (Cd).

Poly- and Perfluoroalkyl Statement:

This product is free from Poly-and Perfluoroalkyl substance (PFAS) such as Perfluorooctanoic and acid (PFOA) and Perfluorooctane sulfonate(PFOS). Raw materials do not contain these substances, and none are intentionally added in manufacturing.