

United Scientific™ SteriTrue™ PS Sampling Scoops

Specifications

Description:

- SteriTrue™ Sterile Sampling Scoop, Polystyrene, 250 mL (8 oz)

Features:


- Ergonomically designed with voice-of-customer input
- Flat bottom surface for stability during weighing
- Gamma sterilized, individually packaged
- Lot stamped for traceability
- PS material meets FDA guidelines

Physical Properties:

- Material: High-impact white polystyrene
- Meets FDA criteria in 21 CFR 177.1520 (2017)
- Heavy Metals free, BPA & Phthalate free, BSE/TSE compliant
- Packaging: sealed polyethylene bag
- Sterilization: Gamma radiation
- Expiration date: 2 yrs from sterilization date
- Overall Length: 233.50mm (9.19 inches)
- Handle Length: 111.50mm (4.39 inches)
- Maximum Width: 72.86mm (2.87 inches)
- Weight: 1.43oz (40.6g)
- Quantity per pack: 100/Pk and 10/Pk
- Country of Origin: India



Item #	Description	Quantity per pack
UNST102	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 250 mL (8 oz)	100
UNST202	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 250 mL (8 oz)	10

 **WARNING:** This product can expose you to chemicals including Styrene, which is [are] known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.