

United Scientific™ SteriTrue™ PS Sampling Scoops

Specifications

Description:

- SteriTrue™ Sterile Sampling Scoop, Polystyrene, 60 mL (2 oz)

Features:


- Ergonomically designed with voice-of-customer input
- Flat bottom surface for stability during weighing
- Gamma sterilized, individually packaged
- Lot stamped for traceability
- PS material meets FDA guidelines
-

Physical Properties:

- Material: High-impact polystyrene
- Meets FDA criteria in 21 CFR 177.1520 (2017)
- Heavy Metals free, BPA & Phthalate free, BSE/TSE compliant
- Packaging: sealed polyethylene bag
- Sterilization: Gamma radiation
- Expiration date: 2 yrs from sterilization date
- Overall Length: 202.20mm (7.96 inches)
- Handle Length: 126.30mm (4.97 inches)
- Maximum Width: 45.00mm (1.77 inches)
- Weight: .78oz (22.2g)
- Quantity per pack: 100/Pk and 10/Pk
- Country of Origin: India



Item #	Description	Quantity per pack
UNST100	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 60 mL (2 oz)	100
UNST200	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 60 mL (2 oz)	10

 **WARNING:** This product can expose you to chemicals including Styrene, which is [are] known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.