

# United Scientific™ SteriTrue™ PS Sampling Scoops



United Scientific™ SteriTrue™ PS sterile sampling scoops are unique designs developed with end-users!

The polystyrene scoops are ergonomically designed for long periods of use, they can hold stated volumes\*, and each scoop has a flat bottom surface for stability during weighing. The most unique feature of the scoops is the long handle: the cross ribs add structural support to help prevent breakage, and the bottom curved design fits comfortably in the palm of the hand.

- **Ergonomically designed with voice-of-customer input**
- **Flat bottom surface for stability during weighing**
- **Gamma sterilized, individually packaged**
- **Lot stamped for traceability**
- **PS material meets FDA guidelines and is safe for food, pharma, and other uses**



**SteriTrue™ Sterile Sampling Scoops are available in three sizes.**

\* Volumes are estimated and dependent on the sampling material. Stated volumes may require a rounded full scoop.

Item #	Description	Quantity per pack
UNST100	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 60 mL (2 oz)	100
UNST200	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 60 mL (2 oz)	10
UNST101	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 125 mL (4 oz)	100
UNST201	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 125 mL (4 oz)	10
UNST102	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 250 mL (8 oz)	100
UNST202	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 250 mL (8 oz)	10

**WARNING:** This product can expose you to chemicals including Styrene, which is [are] known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



**United Scientific®**

1001 Technology Way  
Libertyville, IL 60048

[www.unitedsci.com](http://www.unitedsci.com)  
847-336-7556