

United Scientific™ UniSafe™ Narrow Mouth Vented Wash Bottles, LDPE



These bottles are made of translucent low-density polyethylene. Solutions can be dispensed easily and uniformly by lightly squeezing the bottle. The top closure and stem are molded in one piece to prevent against leaking. The closure is vented which also helps against pressure build-up and leakage.

- Industrial quality low density polyethylene (LDPE) bottle
- Polypropylene (PP) closure
- Leakproof one-piece closure design
- Vented closure to prevent pressure build-up



UN366032
UniSafe™ Narrow Mouth Vented Wash Bottle, LDPE, 500 mL

United Scientific™ UniSafe™ Narrow Mouth Vented Wash Bottles, LDPE

Item #	Capacity (mL)	Height (mm)	Bottle O.D. (mm)	Neck I.D. (mm)	Cap (mm-thread)	Quantity per pack
UN366030	125	151	51	17.5	24 - 415	6
UN366031	250	189	62	17.5	24 - 415	6
UN366032	500	218	73	17.5	24 - 415	6
UN366033	1000	285	92	27	38 - 430	4



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084