

United Scientific™ EcoPoint™ Variable Volume Pipettors



Accurate ergonomic pipetting at a great value!

Our EcoPoint™ variable volume pipettors are designed to dispense liquids with a high degree of accuracy, and are very easy to work with. The lightweight plunger is designed to facilitate use for extended periods of time.

Each pipette box includes a calibrated pipettor, sample pipette tips, a calibration tool, a quality control calibration certificate, and a user's manual.

- Available in these sizes (µL): 0.5-10, 2-20, 5-50, 10-100, 20-200, 100-1000, 1000-5000, 2000-10000
- Designed for use with universal pipette tips (depending on the volume range)
- Fully autoclavable (121°C, 0.15Mpa)
- Manufactured, calibrated, and tested per ISO standards.
- CE certified



EVP-10 EVP-20 EVP-50 EVP-100



EVP-200 EVP-1000 EVP-5000 EVP-10000

United Scientific™ EcoPoint™ Variable Volume Pipettors

Item #	Description
EVP10	United Scientific™ EcoPoint™, 0.5-10 µL
EVP20	United Scientific™ EcoPoint™, 2-20 µL
EVP50	United Scientific™ EcoPoint™, 5-50 µL
EVP100	United Scientific™ EcoPoint™, 10-100 µL
EVP200	United Scientific™ EcoPoint™, 20-200 µL
EVP1000	United Scientific™ EcoPoint™, 100-1000 µL
EVP5000	United Scientific™ EcoPoint™, 1000-5000 µL
EVP10000	United Scientific™ EcoPoint™, 2000-10000 µL



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084

United Scientific™ EcoPoint™ Variable Volume Pipettors



Accurate volume
adjust knob

Plunger with
smooth movement

Ergonomic
finger hook

Quick release
tip ejector

Easy to read
display

Universal fit
tipcone compatible
with most tips



What's in the box:

- Calibrated pipettor
- Sample pipette tips
- Calibration tool
- Quality control calibration certificate
- User's manual



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084