

United Scientific™ Microscope Slide Box, ABS



Designed as a convenient solution for temporary storage and transferring of microscope slides. The bottom is lined with a cork pad, and the lid has a printed index for quick reference.

Accommodates microscope slides in two sizes:
 75.0 mm × 25.0 mm × (0.8-1.2 mm)
 76.0 mm × 26.0 mm × (0.8-1.2 mm)

- **Made of strong ABS material**
- **Snap closures and hinged-lids provide safe operation**
- **A water-proof index card is mounted inside the lids for simple specimen identification**
- **Cork sheet protects the slides**
- **The 100-place slide box includes a metal snap closure**



United Scientific™ Microscope Slide Box, ABS

Item #	Description	Quantity
62607	United Scientific™ 25-Place Microscope Slide Box, ABS, White	1
62608	United Scientific™ 25-Place Microscope Slide Box, ABS, Blue	1
62609	United Scientific™ 50-Place Microscope Slide Box, ABS, White	1
62610	United Scientific™ 50-Place Microscope Slide Box, ABS, Blue	1
62611	United Scientific™ 100-Place Microscope Slide Box, ABS, White	1
62612	United Scientific™ 100-Place Microscope Slide Box, ABS, Blue	1



1001 Technology Way
 Libertyville, IL 60048

www.unitedsci.com
 847-336-7556

novatech
 INTERNATIONAL

sales@novatech-usa.com
 www.novatech-usa.com
 Tel: (866) 433-6682 Fax: (866) 433-6684
 Tel: (281) 359-8538 Fax: (281) 359-0084

United Scientific™ Microscope Slide Box, ABS



25-Place Microscope Slide Box



50-Place Microscope Slide Box



100-Place Microscope Slide Box



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084