

United Scientific™ PP Slide Mailers



A convenient and safe solution for holding and transferring microscope slides. The inside walls feature slotted channels to keep slides separated and to make removal easier. The lids close securely to keep the slides safe.

Accommodates microscope slides in two sizes:
75.0 mm × 25.0 mm × (0.8-1.2 mm)
76.0 mm × 26.0 mm × (0.8-1.2 mm)

- **Securely hold slides: 3-place option with screw cap, and 5-place option with hinge cap**
- **Durable Polypropylene**
- **Holds two common slide sizes**

United Scientific™ PP Slide Mailers

Item #	Description	Quantity per Pack
62605	United Scientific™ 3-Place Slide Mailer with Screw Cap	100
62606	United Scientific™ 5-Place Slide Mailer with Hinge Cap	100



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084