

## Durable & Reusable

- Constructed of durable ABS material
- Easy to clean and re-use for multiple projects

## Secure & Safe

- Grooved Slots to securely hold slide in place
- Lining material for added cushion against breakage
- Locking feature on outside of box to prevent accidental opening

## Systematic

- Multiple identifying features to make labeling and classification simple
- Removable slide inventory sheet included for slide identification

## Microscope Slide Box

- Lockable clasp and hinge pins are nickel-plated and rust resistant
- Identification made easy with numbered cork lining corresponding directly to the slide sheet. Outside cover also features an extra labelling area for "Project" and "Date" classification
- Stack multiple boxes with molded in feet to steadily hold for transport or storage



Cat. No.	Capacity	Color	Material	Lining	Dimensions (L x W x H)	
					Inches	Centimeters
1400-SLCS201	100 slides	Blue	ABS	Cork	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4

## Premium Microscope Slide Box Assortment

- All plastic construction eliminates the potential of rust
- Sliding polypropylene lock has extra security feature allowing use of a 2-mm zip tie or security tag (not included)
- Center bar in slide storage area is clearly numbered for easy correspondence to the slide inventory sheet. Assorted colors can be used for additional project coding
- Sold as a pack of five slide boxes with assorted colors: one each of Blue, Orange, Green, Red, and Purple



Cat. No.	Capacity	Color	Material	Lining	Dimensions (L x W x H)	
					Inches	Centimeters
1400-SLCS101	100 slides, each box	Blue, Orange, Green, Red, and Purple	Translucent ABS	EVA Foam	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2