







LABO50LCD

DATA SHEET

Head	6-inch LCD screen (640 x 480 pixels resolu-
	tion) with 2.0 MP digital camera, with SD card
	slot, USB input and AV port.
Nosepiece	Quadruple
Eyepieces	4X, 10X, 40X (S) and 100X (S) (oil), metallic
	finish, anti-slip ring and colored ring for easy
	identificationN/A
Optical System	Achromatic lenses
Stage	Double layer with X-Y coaxial movements,
	110 X 110 mm.
Condenser	Abbe, A. N. 1.25
Diaphragm	Iris type with filter holder.
Focus	Coarse an d f ine anti-slip control knob
llumination	LED with intensity adjustment.
Dimensions	***
Power Supply	110-240V, 50/60 Hz AC adapter.
Warranty	5 years.
•	

LABOTECA BY VELAB, Co. offers Limited Warranty (Parts and Labor) for components that fail due to defects in materials or assembly. The guarantee starts from the delivery date.

The warrantee will be canceled if the equipment is not returned in the original packaging and the correct processed claim documentation. All claims are the sole discretion of VELAB, Co.

VELAB, Co. reserves the right to make changes in the technology, features, specifications and design of its equipment without prior notice.