

## LABO Z1

### DATA SHEET



<b>Ocular</b>	WF 10X /18 mm with block screw.
<b>Nosepiece</b>	N/A
<b>Eyepieces</b>	2X (20 total increases)
<b>Optical System</b>	N/A
<b>Stage</b>	Interchangeable, double-sided white / black with a pair of tweezers.
<b>Condenser</b>	N/A
<b>Diaphragm</b>	N/A
<b>Focus</b>	Coarse anti-slip control knobs.
<b>Illumination</b>	LED.
<b>Dimensions</b>	120 x 250 x 130 mm / 4.7" x 9.8" x 5.1"
<b>Power Supply</b>	Two 1. 2V AAA batteries.
<b>Warranty</b>	5 years.

LABOTECA BY VELAB, Co. offers Limited Warranty (Parts and Labor) for components that fail due to defects in materials or assembly. The guarantee starts from the delivery date.

The warranty will be canceled if the equipment is not returned in the original packaging and the correct processed claim documentation. All claims are the sole discretion of VELAB, Co.

**VELAB, Co. reserves the right to make changes in the technology, features, specifications and design of its equipment without prior notice.**