

LABOTECA



VELAB Co.



1431 West Polk Ave.
Pharr, TX 78577, USA



+1 956 429 3410



sales@velab.net



www.velab.net

WORLD



LABOTECA


BINOCULAR BIOLOGICAL MICROSCOPE

Used in the educational field, recommended for middle school and high school biology classes. Its LED illumination saves energy, has a long lifespan and allows comfortable observation. It is a portable device that can be used in laboratories or outdoor operations using the rechargeable batteries included

A: Includes rechargeable batteries, dust cover, Instruction manual and power supply adapter

B: 10-year warranty.

C: Practical and easy to operate design. Used for dissection of plants, insects, fungi and basic observation of cells.



TECHNICAL DATA

Head	Monocular head 45° inclined, 360° rotatable.
Eyepieces	Pair WF 10x -18 mm.
Nosepieces	Triple revolving nosepiece, 360° full rotation..
Objectives	Achromatic Objectives 4X, 10X, 40X. With non-skid color ring for easy identification.
Body	Stationary, robust and reinforced body.
Stage	Flat stage with two clips, 110 x 110 mm..
Condenser	Abbe A.N. 0.65. LED
Diaphragm	Iris diaphragm.
Focusing	Coaxial coarse and calibrated fine focus.
Illumination	LED lamp with light intensity adjustment
Power	Power supply 120 V-50/60 Hz, includes 3-AA rechargeable batteries.