PB SERIES

Readability 0.01 g





VE-1000

Features

- · High contrast LCD screen.
- · Large stainless steel pan.



















- - · Capacity tracker
- · Leveling feet: Two back feet
- · RS 232 interface.
- · Multiple weighing units.

Standard Functions

- · Weighing
- · Selectable weighing
- · Percentage of weight
- · Piece counting"

Accesories Included

- · Calibration weight 100g and 200g · RS232 cable
- · Transparent protective cover
- · Dust cover

- · USB adapter
- · Hook for suspended weighing"

TECHNICAL DATA

Model **VE-1000** Capacity 1000 g Readability 0.01 g Repeatability ±0.01 g Linearity ±0.02 g Weighing Units g, ct, oz, ozt, dwt, GN, lb, N, dr, tlT, tlJ, tlH, T, T/A/R, /A/R, kg Calibration External. Pan Size Ø 133mm / Ø5.2" Stabilization Time ≤ 3 seconds Interface RS-232, only can comunicate with computer Display size 117 X 35mm, digit height 20mm / 0.8" Housing ABS plastic. **Operating Environment** 5° to 35°C / 41° to 35°F / Humidity 20% - 80% **Overall Dimensions** 295 x 208 x 88 mm / 11.6" x 8.1" x 3.4" Underhook mechanism Hook included **Power Supply** 110/240V, 50/60 Hz AC adapter and rechargeable battery. Net Weight 3.6 kg Warranty 1 year.



www.velab.net

+1 956.429.3410

► sales@velab.net

▶ 1431 W Polk Ave, Pharr Texas