PRECISION BALANCES

H SERIES

GLP GUIDELINES

Readability **0.1 g**



VE-15001H VE-20001H

VE-30001H

Features

- High contrast LCD screen.
- Robust ABS housing.
- Leveling feet: Two back feet.
- Multiple weighing units (24 total).

Standard Functions

- Weighing.
- Parts counting.
- Percentage weighing.

Accesories Included

• Rechargeable battery.

TECHNICAL DATA

Model	VE-15001H	VE-20001H	VE-30001H
Capacity	15,000 g	20,000 g	30,000 g
Readability	0.1 g	0.1 g	0.1 g
Repeatability	±0.2 g	±0.2 g	±0.2 g
Linearity	±0.2 g.	±0.2 g.	±0.2 g.
Weighing Units	g, ct, lb, oz, dr, GN, twt, kg (24 in total).		
Calibration	External.		
Pan Size	330 x 260 mm / 13" x 10"		
Stabilization Time	2 seconds.		
nterface	Bluetooth		
Display	LCD, digit height 25mm / 1"		
Housing	ABS plastic.		
Operating Environment	10° to 30°C / 50° to 85°F / Humidity 20% - 80%		
Overall Dimensions	293 x 219 x 80.5 mm / 155" x 13" x 3.5		
Power Supply	CA 110/240V, 50/60 Hz and rechargeable battery.		
Net Weight	7.5kg / 17 lb.		
Warranty	l year.		



www.velab.net

- +1 956.429.3410
- ▶ sales@velab.net
- ▶ 1431 W Polk Ave,Pharr Texas



- Internal rechargeable battery.
- \cdot Large stainless steel pan.
- \cdot Bluetooth