## **VE-153G**

PRIME

Zoom System, 0.7X, 1X, 1.5X, 2X, 3X, 4X y 5X



### Main Features

- · Head with mini USB, USB, HDMI ports.
- · Adjustable screen -5 ° to 15 ° inclined.
- · Stative, robust and reinforced.
- · SD card slot.
- · Built-in SOPTOP software for snapshot editing.









- · Software Velab, helpful for digital editing and measure in realtime, setting image parameters and general features.
- · HotSpotViewWiFi app, compatible with Android and iOS, for viewing images y capturing instantaneously on other devices.

# **Applications**

- · Education: Upper grade levels
- · Industrial

· Research

## Accesories Included

· USB cable

· HDMI cable

· Power cord

· Software

#### TECHNICAL DATA

**VE-153G** Model

Fixed screen 10", IPS with HD resolution, 4.0 MP. Camera

0.7X, 1X, 1.5X, 2X, 3X, 4X and 5X. Objectives

Zoom system. Optical System

Fixed. Stage

Coarse anti-slip. Focus

Independent LED ring, temperature 5000-5500K. Ilumination

320 x 482 x 260 mm / 10.2" x 12.5" x 18.2". Dimensions

110/240V, 50/60 Hz AC adapter. **Power Supply** 

10 years Warranty

VELAR

www.velab.net

**+**1 956.429.3410

▶ sales@velab.net

▶ 1431 W Polk Ave, Pharr Texas 78577