



VE-EL



► **Features**

- Easy to use.
- Robust construction for demanding use
- Retractable height measurement bar

► **Standard Functions**

- Weighing

► **Accessories Included**

- n/a

► **TECHNICAL DATA**

	EL
Model	
Capacity	200 kg /440 lb
Readability	0.1kg / 0.2lb
Repeatability	±0.2kg / 0.4lb
Linearity	±0.4lb / 0.2kg
Weighing Units	g, kg, lb
Calibration	External
Pan Size	405 X 280 mm. / 15.94" x 11"
Stabilization Time	4 seconds
Housing	Steel
Operating Environment	10° to 30°C / 50° to 85°F / Humidity 20% - 80%
Overall Dimensions	630 x 465 x 148 mm / 24.80" x 18.30" x 5.82"
Net Weight	25.5 kg / 56.21 lbs
Warranty	1 year