



► Features

- High performance camera
- Software to capture ImageView, analyze and edit captures
- Capable of making measurements in a wide variety of units: pixels, nanometers, micrometers, millimeters, centimeters, meters, inches.

► Applications

- Research laboratories
- Education: Upper level

► Accessories included

- 2 adapter rings 30.5 mm
- USB cable
- Adapter for microscope FMA050 0.5X
- Software

TECHNICAL DATA

Model	VE-LX1800
Use	Video camera for PC.
Total pixels	18 MP
Image sensor	1/2.3" CMOS
Resolution	4912 x 3684 MP, 2456 x 1842 MP, 1228 x 922 MP.
White balance	ROI white balance / manual tint temperature adjustment / NA for the monochrome sensor
Interface	USB 3.0
Operating system	Microsoft Windows XP, Vista, 7/8./10 (32 & 64 bits) (Mac OS X) and Linux.
Recording system	Imagen and video.
Power supply	DC 5V via USB cable
Warranty	1 year