



► Features

- Real-time image observation on computer
- Capture of images for storage

► Applications

- Research laboratories
- Education.

► Accesories included

- Adapter rings 30 and 30.5 mm
- Software
- USB cable

TECHNICAL DATA

Model	VE-MC5
Use	Video camera for PC.
Total pixels	5 MP
Image sensor	1/2" CMOS
Pixel size	3.2 μm x 3.2 μm
Video output	24 bit color
Resolution	640 x 480, 1024 x 768, 2048 x 1536, 2530 x 1978 MP
White balance	Automatic and manual
Brightness	1.8 U / LUX-SEG 550 NM
Interface	USB 2.0
Power supply	CD 5V vía cable USB
Warranty	1 año