



► Main Features

- Binocular head, inclined at 60 °.
- Interpupillary distance of 55 ~ 75mm
- Stable, robust, reinforced and easy to transport.
- Lighting with independent switches.
- Diopter adjustment in one eyepiece.
- Anti-slip ring on lens and cap.
- Stage with clamps.

► Applications

- Education: Lower grade level.
- Research
- Industrial

► Accessories Included

- Power cord
- Dust cover

TECHNICAL DATA

Modelo	VE-S4
Ocular	WF 10X/20 mm.
Objectives	Variable Magnification Zoom System 1X, 1.5X, 2X, 2.5X, 3X, 3.5X, 4.5X (10 to 45 total increases).
Stage	Interchangeable, translucent white and white / black.
Focus	Coarse anti-slip control knob.
Illumination	Transmitted and incident LED.
Dimensions	200 x 340 x 275 mm / 7.8" x 13.3" x 10.8".
Power Supply	110/240V, 50/60 Hz AC adapter.
Warranty	10 years