



► Main Features



- Diopter adjustment on both eyepieces.
- Binocular head with interpupillary distance adjustment.
- Head inclined at 45°.
- Stage with clips.
- Software for camera operating, capturing photo, video and digital analysis.
- Stative, robust, reinforced and easy to transport.
- Zoom adjustment knobs.
- Interpupillary distance of 55 ~ 75 mm.
- Integrated digital camera.
- Lighting with independent intensity adjustment.

► Applications

- Education: Upper grade levels
- Industrial
- Research

► Accessories Included

- USB cable
- Software
- Dust cover

TECHNICAL DATA

Model	VE-S5C
Eyepiece	WF 10X/20 mm.
Camera	1.3 MP, integrated.
Objectives	0.7X, 0.8X, 1X, 1.5X, 2X, 3X, 4X, 4.5X (7 to 45 total increases).
Optical System	Greenough System, with variable magnification zoom.
Stage	Interchangeable, translucent white and white / black.
Focus	Macrométrico antiderrapante.
Illumination	Transmitted and incident LED.
Dimensions	225 x 355 x 270 mm/ 8.8" x 13.9" x 10.6".
Power Supply	CA 110/240, 50/60 Hz.
Warranty	10 years