



► Main Features

- Binocular head, inclined at 45 °.
- Interpupillary distance of 55 ~ 75mm
- Stable, robust, reinforced and easy to transport.
- Lighting with independent switches.
- Diopter adjustment in both eyepieces.
- Anti-slip ring on lens and cap.
- Stage with clamps.

► Applications

- Education: Lower grade level
- Research
- Industrial

► Accesories Included

- Power cord
- Dust cover

► TECHNICAL DATA

Modelo	VE-S6
Ocular	WF 10X/22 mm.
Objectives	Variable Magnification Zoom System 0.67X, 0.7X, 0.8X, 1X, 1.5X, 2X, 3X, 4X y 4.5X (6.7 to 45 total increases).
Stage	Interchangeable, translucent white and white / black.
Focus	Coarse anti-slip control knob.
Illumination	Transmitted and incident LED.
Dimensions	245 x 435 x 312 mm/ 9.6" x 17.1" x 12.2".
Power Supply	110/240V, 50/60 Hz AC adapter.
Warranty	10 years