

## X-ZOOMI MONOCULAR MICROSCOPE

The X-ZOOMI monocular microscope was designed to introduce children to the world of science. Children will spend hours of fun learning how to observe all sorts of samples, making the learning process an enjoyable adventure into the world of microscopy.



### CHARACTERISTICS

**EYEPIECES:** WF 10X / 18mm with stud

**HEAD:** Monocular

**NOSEPIECE:** triple nosepiece with click-stops

**OBJECTIVES:** achromatic 4X, 10X and 40X, with color ring for easy identification.

**STAGE:** 90 x 90 mm, with 2 clamps.

**ILLUMINATION:** LED transmitted (base) and incident (arm)

**CONDENSER:** NA 0.65

**DIAPHRAGM:** Disc

**DIMENSIONS:** 11.43 x 16.51 x 29.21 cm

**POWER SUPPLY:** 3 AA 1.5V batteries (not included).

---

### ACCESSORIES:

- 5 fixed preparations
- 1 slide with 1 cavity
- Slide labels
- Coverslips
- Plastic dropper with red and blue dye
- Tweezers
- Tissue paper
- Petri dish
- Scalpel
- Dissection needle
- 1 transfer pipette
- 1 test tube
- Plastic sheath