


# VEE GEE®

VEE GEE Model 3020 Double Buret Clamps hold two burets from micro to 100ml size. Constructed of chrome plated, zinc-aluminum alloy for durability, the unit is simple to operate -- compress the scissor-like mechanism, insert buret, and release. Numbers and graduations on burets are easy to read. Attach easily to support rods with a convenient hook clamp. Knob for tightening is on the front for easy access.



## Clamps

VEE GEE®, Buret, Double

Cat. No.	Style	Max. Buret Size
3020	Double	100ml



### VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
 www.veegee.com • email: sales@veegee.com

**novatech**  
 INTERNATIONAL

sales@novatech-usa.com  
 www.novatech-usa.com  
 Tel: (866) 433-6682 Fax: (866) 433-6684  
 Tel: (281) 359-8538 Fax: (281) 359-0084

No.VGBRC081913/3020  
 Printed in the USA