



VEE GEE Digital Thermometer Series:

- High Accuracy Standard Model
- Ultra Accurate Model
- 12" Long Probe Model

VEE GEE® digital stem thermometers are ideal for temperature monitoring in most laboratory settings. Great for monitoring reagents and liquids in test tubes, beakers, & flasks, VEE GEE digital thermometers are an essential tool used in food processing, refrigeration, and horticulture.

High Accuracy, Cat No's: **83210**, (127 mm / 5" probe) & **83210-12**, (302mm / 12" probe)

Range: -50 to +150°C / -58 to +302°F

Resolution / Accuracy: 0.1°/±1.0°C & 0.1°/±2.0°F

The Ultra Accurate, Cat No. **83110**, (127 mm / 5" probe)

Range: 10 to +100°C / +14 to +212°F

Resolution / Accuracy: 0.1°/±0.5°C & ±1.0°F

The bi-metal stem is 3.5mm (0.138") in diameter and includes:

- a protective case featuring a pocket clip
- a 1.5V alkaline coin cell battery (Type LR44)
- automatic power-off after 10 minutes of idle use
- 180° swivel head
- Face Dimensions 36 X 17mm/1.4 X 0.7 in.



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegoe.com • email: sales@veegoe.com

VBRC20216/83000-Digital-Series
Printed in the USA