

Versilon™ C-210-A Tubing

Oil, Fuel and Grease Resistant

C-210-A Tubing

The Versilon C-210-A tubing provides exceptional stability when exposed to oils, fuels, and grease. Use this tough ether-based polyurethane tubing for transferring abrasive and viscous slurries, dispensing lubricants and degreasers, and recovering coolants. The plasticizer-free tubing offers excellent chemical resistance and withstands temperatures ranging from -73 to 79 °C (-100 to 175 °F).



CERTIFICATION

Cole-Parmer will provide a certificate of compliance FREE at your request where applicable.

APPLICATIONS

- General-purpose, laboratory, and single-use applications

FEATURES/BENEFITS

- Free technical application support available to assist with tubing selection
- Excellent abrasion and tear resistance
- Remains flexible even in subzero temperatures
- Consistently tight lot-to-lot tolerances
- Clear tubing to view fluid being transferred

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
France	+33 (0) 1486 37800

Germany	+49 (0) 9377 92030
India	+1.800.266.1244
Italy	+39 (0) 1313 89513
UK	+44 (0) 1480.272279
All others	+1.847.549.7600

Specification and Ordering Table

Dimensions: in. (mm)			Max pressure: PSI (bar)	Item number	Roll length (ft)
Tubing ID	Tubing OD	Wall			
7/16"	11/16"	0.125 (3.2)	49 (3.4)	50107-04	50
1/4"	5/16"	0.03125 (0.8)	28 (1.9)	50106-95	50
1/2"	3/4"	0.125 (3.15)	46 (3.2)	50107-07	50
5/8"	3/4"	0.0625 (1.55)	24 (1.6)	50107-08	50
5/8"	13/16"	0.09375 (2.35)	32 (2.2)	50107-09	50
5/16"	7/16"	0.0625 (1.55)	36 (2.45)	50106-99	50
1/4"	1/2"	0.125 (3.15)	70 (4.8)	50106-98	50
7/16"	5/8"	0.09375 (2.4)	40 (2.7)	50107-03	50
3/8"	1/2"	0.0625 (1.6)	34 (2.3)	50107-00	50
1/4"	7/16"	0.09375 (2.35)	58 (4)	50106-97	50
3/16"	5/16"	0.0625 (1.6)	56 (3.8)	50106-93	50
3/16"	3/8"	0.09375 (2.35)	70 (4.8)	50106-94	50
1/8"	1/4"	0.0625 (1.6)	74 (5.1)	50106-91	50
3/16"	1/4"	0.03125 (0.8)	34 (2.3)	50106-92	50
3/8"	5/8"	0.125 (3.2)	54 (3.7)	50107-02	50
1/16"	1/8"	0.03125 (0.8)	70 (4.8)	50106-89	50
3/8"	9/16"	0.09375 (2.4)	45 (3.1)	50107-01	50
1/2"	5/8"	0.0625 (1.6)	26 (1.8)	50107-05	50
1/8"	3/16"	0.03125 (0.8)	45 (3.1)	50106-90	50
1/4"	3/8"	0.0625 (1.55)	42 (2.9)	50106-96	50
1/2"	11/16"	0.09375 (2.4)	36 (2.5)	50107-06	50
5/8"	7/8"	0.125 (3.15)	38 (2.6)	50107-10	50
3/4"	15/16"	0.09375 (2.4)	26 (1.8)	50107-11	50
1"	1-1/4"	0.125 (3.2)	28 (1.9)	50107-13	50
3/4"	1"	0.125 (3.2)	33 (2.3)	50107-12	50



sales@novatech-usa.com

www.novatech-usa.com

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084

Cole-Parmer®
an antylia scientific company